

VICAIR MATTRESS 415



SKU: ARM

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Pain relief, unrest reduction, pressure ulcer prevention. These are just some examples of the immediate relief provided by the Vicair 415 Mattress.

The design offers optimal pressure redistribution, reduction of shear forces, optimal immersion, and stable, yet comfortable, positioning. Even with existing decubitus stage I-IV (clinically proven). If a situation requests an individual solution: custom sizes and shapes are available on request.

The Vicair SmartCell air filling makes the mattress an extremely lightweight and non powered alternative to other air technique replacement mattresses.

FEATURES

- Optimal pressure redistribution
- Pain relief
- Extremely suitable for immobile users

SPECIFICATIONS

Height	15 cm
Average weight	10 kg
Max. user weight	250 kg
Ready for use	Yes
Standard sizes	3 (Size Chart)
Machine washable	No

Latex free
Warranty

Yes
5 years

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au