



VICAIR MULTIFUNCTIONAL O2 WHEELCHAIR BACK CUSHION

\$1,169.10 - \$1,751.00 Price range: \$1,169.10 through \$1,751.00



SKU: OMFB

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Excellent pressure redistribution, comfort, and support make the Vicair Multifunctional O2 back support cushion perfect for use in a large variety of situations. It can be used in several types of back systems and in tilt-in-space chairs with adjustable seat depth. The cushion is easily fixed onto the solid back support by use of hook and loop strips.

With the Vicair Back Connector, the Multifunctional O2 Back can also be used in tension adjustable backs such as the Vicair Strap Back. Although it automatically adapts to the body's contours, the cushion can also be adjusted for optimal positioning, while the user is still leaning against the cushion. 100% machine washable and breathable, for optimal hygiene and perfect temperature and moisture regulation.

The Vicair Multifunctional O2 back helps prevent pressure ulcers. It is also a high quality addition to a pressure ulcer treatment plan (stage I-IV).

FEATURES

- Optimal pressure redistribution helps prevent pressure ulcers
- Top compartments support the shoulders | compartment shape provides freedom of movement
- Multidimensional lumbar support
- Bottom compartment adjustable for sacral support
- Multidimensional thoracic support
- Supporting the spine offloading if required
- Fully machine washable
- SWL: 150kg

SPECIFICATIONS

Length	36-60 cm
Standard sizes	7 (Size Chart)
Average weight	975 grams
Loss of seat depth	3 cm
Max. user weight	150 kg
Machine washable	Yes. Max. 60C
Latex free	Yes
Warranty	24 months on cushion's inner cover Lifetime Warranty on SmartCells