

VICAIR O2 BARIATRIC CUSHION WITH COVER

\$1,176.00 - \$1,949.00 Price range: \$1,176.00 through \$1,949.00



SKU: OBA/CI/5555-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **Vicair O2 Bariatric Cushion with Cover** is specifically developed for bariatric wheelchair users requiring advanced pressure care and comfort. Designed to support individuals weighing up to 500 kg, it provides exceptional skin protection, stability, and postural positioning while maintaining optimal hygiene and temperature control.

Features

- Developed specifically for bariatric users, supporting up to 500 kg
- Filled with larger SmartCells for improved stability and pressure distribution
- Based on the Vicair Adjuster O2 design, offering excellent positioning and comfort
- Automatically manages mild pelvic obliquities, with adjustable support for more extreme asymmetries
- 100% machine washable and breathable for superior hygiene and moisture regulation
- Extremely lightweight despite its large size
- Available with multiple cover options including Incotec, Comfair, and Top Cover (2 Pack)

Specifications

Size (mm)	Cover Type
550 x 550	Incotec Cover
550 x 500	Top Cover (2 Pack)
750 x 500	Comfair Cover
750 x 500	Incotec Cover
750 x 500	Top Cover (2 Pack)
650 x 500	Top Cover (2 Pack)
710 x 560	Incotec Cover
800 x 510	Top Cover