

VICAIR POMMEL O2 WHEELCHAIR CUSHION





SKU: OPC10

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Vicair Pommel O2 cushion is ideal for individuals with adducted conditions, where the knees or thighs are pulled together. The cushion's detachable medial thigh support gently separates the knees and/or thighs relieving the area of pressure that could normally result in pain. The SmartCell filling ensures a softer padded surface between the knees and/or thighs, redistributing pressure in this area and therefore also reducing the risk for pressure ulcers. The detachable medial thigh support has its own separate cover, so it can be removed to facilitate forward transfers or catheterisation whilst remaining in the wheelchair.

FEATURES

- Optimal temperature and moisture regulation
- Fully machine washable
- Ideal for users with adducted conditions
- Detachable medial thigh support
- Unique: cushion and medial thigh support covered separately for ease of use

COVER TYPES

Vicair Top Cushion Cover (Qty 2)

- Ideal to help manage one of the main causes of pressure sores
- Standard cover for the Vicair O2 cushions
- Open bottom allows the cushion to breathe continuously
- Optimal ventilation and comfort



Machine washable

Vicair Comfair Cushion Cover

- Made out of multi-dimensional stretch material
- Ergonomically shaped cover for better envelopment and immersion
- Air transmission and heat dissipation
- Machine washable

Vicair Incotec Cushion Cover

- Usable in all situations where incontinence and hygiene is an issue
- Fluid and incontinence resistant
- Extremely usable for institutional and home use
- Machine washable

SPECIFICATIONS

Available height 10 cm

Average weight 850 grams - depending on size

Max. user weight 250 kg Ready for use Yes

Standard sizes 65 (<u>Size Chart</u>)

Machine washable Yes. Max. 60°C

Latex free Yes

Warranty 24 months

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

revnella@medimart.com.au



VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au