


VICAIR POMMEL O2 WHEELCHAIR CUSHION

\$1,258.00 - \$1,373.00 Price range: \$1,258.00 through \$1,373.00



SKU: OPC10

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Vicair Pommel O2 cushion is ideal for individuals with adducted conditions, where the knees or thighs are pulled together. The cushion's detachable medial thigh support gently separates the knees and/or thighs relieving the area of pressure that could normally result in pain. The SmartCell filling ensures a softer padded surface between the knees and/or thighs, redistributing pressure in this area and therefore also reducing the risk for pressure ulcers. The detachable medial thigh support has its own separate cover, so it can be removed to facilitate forward transfers or catheterisation whilst remaining in the wheelchair.

FEATURES

- Optimal temperature and moisture regulation
- Fully machine washable
- Ideal for users with adducted conditions
- Detachable medial thigh support
- Unique: cushion and medial thigh support covered separately for ease of use

COVER TYPES

Vicair Top Cushion Cover (Qty 2)

- Ideal to help manage one of the main causes of pressure sores
- Standard cover for the Vicair O2 cushions
- Open bottom allows the cushion to breathe continuously
- Optimal ventilation and comfort

- Machine washable

Vicair Comfair Cushion Cover

- Made out of multi-dimensional stretch material
- Ergonomically shaped cover for better envelopment and immersion
- Air transmission and heat dissipation
- Machine washable

Vicair Incotec Cushion Cover

- Usable in all situations where incontinence and hygiene is an issue
- Fluid and incontinence resistant
- Extremely usable for institutional and home use
- Machine washable

SPECIFICATIONS

Available height	10 cm
Average weight	850 grams – depending on size
Max. user weight	250 kg
Ready for use	Yes
Standard sizes	65 (Size Chart)
Machine washable	Yes. Max. 60°C
Latex free	Yes
Warranty	24 months