



VULKAN HINGED KNEE SUPPORT

\$130.00 Original price was: \$130.00.\$71.00 Current price is: \$71.00.While Stocks LastOnline Only



SKU: VLK-1032-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **Vulkan Hinged Knee Support** delivers targeted collateral stabilisation for mild to moderate knee instabilities, including meniscus injuries, ligament strains or post-trauma rehabilitation. Equipped with discreet hinges and a soft-cushioned felt lining, it supports the knee's natural movement while preventing awkward valgus/varus deviation and hyperextension. Designed for stability and warmth, it also helps improve circulation and proprioception—all with comfort in daily wear.

Features

- Hidden hinges offer lateral support to control knee deviation and prevent hyperextension
- Soft-cushioned felt contours to the knee for extra comfort
- Thermal effect boosts blood circulation and supports healing
- Latex-free 75% neoprene blend with 25% nylon for durability and heat retention
- Easy to use: pull on, align hinges with knee joint, and secure straps for fit

Specifications



• Material: 75% neoprene blend, 25% nylon

• Sizes Available: Small (30–33 cm), Medium (33–36 cm)

• Care Instructions: Wash in cold water, gentle cycle with mild soap; air-dry only. Do not bleach, iron, or expose to heat or tumble-dry.

• Warranty: 3 months