




VULKAN POSTURE BRACE



SKU: VLK-2075-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The VULKAN Posture Brace is designed to provide customised support to the spine in order to promote correct posture. Equipped with aluminium stays located on either side of the spine and made from breathable elastic materials this product is a user-friendly pocket design to ensure easy application.

The VULKAN Posture Brace is ideal for spondylosis, osteoporosis, and kyphosis.

This product is available in two lengths - For correct length measure the height of the back from C7 to sacrum.

FEATURES

- Easy application.
- Bilateral shoulder straps with pads.
- Made of stretch elastic material for flexible compression.

SIZING

Measure around scapula under arm.

XS	46 - 61cm
S	61 - 76cm
M	76 - 91cm
L	91 - 107cm
XL	107 - 122cm
XXL	122 - 137cm
XXXL	137 - 152cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au