



VULKAN SACRO BELT



SKU: VLK-2262-H

PRODUCT DESCRIPTION

FEATURES

- May be used post-partum and following traumatic injury.
- Relieves instability and stress through compression of the sacroiliac joints.
- Reinforcement strap provides additional compression.
- Hook and loop closure for easy application.

HOW TO WEAR

- For use as Maternity Belt Position belt just below belly for support. Bring straps to back and fasten securely. Pull side straps to adjust tightness then fasten.
- For use as Waist Belt 1.Wrap belt around the waist. 2.Bring straps to the front and fasten. 3.Pull side straps and secure to desired compression.

APPLICATION

Pregnancy

SIZING CHART

Measure around waist.



S 86 - 92cm M 92 - 97cm L 97 - 102cm XL 102 - 110cm XXL 110 - 118cm

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468

support@medimarnt.com.au