



DAYS MEAL TRAY WALKER WITH HANDBRAKES

\$253.00



SKU: HLC-109NS-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Days Adjustable Height Trolley Walker is a sturdy and robust walking aid for use within the home, providing a means of transporting items from one room to another safely.

FEATURES

- Includes clip-on plastic trays that can be removed easily for cleaning.
- Can be used with or without trays or it can be used as a table.
- Rails on tray edges prevent items from slipping off.
- Lockable hand brakes mean this trolley can also be used as a walking aid and provides additional safety when standing from a seated position.
- Handles are adjustable in height.
- Easily accessible loop lockable brakes.
- Supplied with a cream painted steel frame.
- Comes in white or beige colour.

SPECIFICATIONS

- Overall Width: 585mm
- Overall Depth: 560mm
- Handle Height Range: 755mm - 905mm
- Wheel Diameter: 150mm (6")
- Maximum User Weight: 150kg