




REDGUM ALUMINIUM KNEE WALKER WITH HANDBRAKES

\$405.00



SKU: RG30KW

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

A knee scooter or knee walker is a four-wheeled alternative to crutches or a traditional walker as a mobility aid.

Today's version is a lightweight, foldable design that, with the knee flexed, supports the shin of the unusable limb.

The knee walker is a safe, comfortable, and easy-to-maneuver alternative to the traditional crutch.

The knee walker is foldable and suitable for transport in the boot of a car.

KEY FEATURES

- Aluminium powder-coated frame.
- Compact & Easy to fold
- Multi height adjustable platform
- Adjustable handlebar height
- Quick-adjust for left or right leg
- Dual locking brakes
- 8" wheels / Carry basket

SPECIFICATIONS

Maximum User Weight:	130kg
Knee Pad Height (from ground level):	40cm - 51.5cm
Handle Height (from ground level):	79cm - 95cm

Product Weight:

12kg

Castor (Wheel) Diameter:

20cm (8")

Buy this Redgum Aluminium Knee Walker online at MediMart using Afterpay or [NDIS method](#) today!