



# REDGUM ALUMINIUM KNEE WALKER WITH HANDBRAKES

\$405.00 Original price was: \$405.00.<u>\$365.00</u>Current price is: \$365.00.



SKU: RG30KW

## **PRODUCT DESCRIPTION**

#### **PURPOSE & BENEFITS**

A knee scooter or knee walker is a four-wheeled alternative to crutches or a traditional walker as a mobility aid.

Today's version is a lightweight, foldable design that, with the knee flexed, supports the shin of the unusable limb.

The knee walker is a safe, comfortable, and easy-to-manoeuvre alternative to the traditional crutch.

The knee walker is foldable and suitable for transport in the boot of a car.

# **KEY FEATURES**

- Aluminium powder-coated frame.
- Compact & Easy to fold
- Multi height adjustable platform
- Adjustable handlebar height
- Quick-adjust for left or right leg
- Dual locking brakes
- 8" wheels / Carry basket

### **SPECIFICATIONS**

Maximum User Weight: 130kg

Knee Pad Height (from ground level): 40cm - 51.5cm

Handle Height (from ground level): 79cm - 95cm



Product Weight: 12kg
Castor (Wheel) Diameter: 20cm (8")

Buy this Redgum Aluminium Knee Walker online at MediMart using Afterpay or NDIS method today!