


## WHILL MODEL C2 POWER WHEELCHAIR

~~\$10,788.00~~ Original price was:

~~\$10,788.00~~ ~~\$9,709.20~~ Current price is: \$9,709.20. Sale

Available Online Only!



SKU: PWS655410-H



## PRODUCT DESCRIPTION

### Purpose & Benefits of the Whill Model C2 Power Wheelchair

The Whill Model C2 Power Wheelchair is a cutting-edge mobility solution, offering advanced manoeuvrability and comfort in a sleek, compact design. Ideal for both indoor and outdoor use, the C2 combines powerful performance with intuitive controls, allowing users to navigate through tight spaces and over various terrains with ease. Its lightweight frame and easy disassembly make it perfect for travel, while the ergonomic seating ensures comfort for extended use. The Whill Model C2 is designed for those seeking independence and freedom, providing a reliable and stylish way to stay mobile.

### Key Features

- **Compact and Lightweight Design:** The Model C2 is designed to be easily disassembled into three lightweight parts, making it ideal for transport and storage.
- **Advanced Driving Control:** Features an intuitive joystick for precise control, allowing users to navigate smoothly through tight spaces and around obstacles.
- **All-Terrain Capability:** Equipped with powerful motors and large front omni-wheels, the C2 can handle various surfaces, including gravel, grass, and uneven pathways, offering reliable outdoor performance.
- **Customisable Seating:** Ergonomic seat with adjustable backrest and footrest for personalised comfort, ensuring optimal support during long periods of use.
- **Remote Control Functionality:** Includes a remote control feature, allowing users to move the wheelchair via smartphone app for easy access in and out of tight spaces.

## Specifications

- **Top Speed:** 10 km/h
- **Range:** Up to 18 km on a single charge
- **Weight Capacity:** 136 kg
- **Turning Radius:** 76 cm, ideal for navigating indoors
- **Weight:** 52.9 kg (without battery)
- **Battery Type:** Lithium-ion, easy to remove and charge
- **Dimensions:** 98 cm (L) x 55 cm (W)
- **Cleaning:** Wipe clean with a damp cloth