



WINDSOR PORTABLE TOILET/COMMODE LIFT SEAT

\$3,501.00 Original price was:

\$3,501.00.<u>\$3,151.00</u>Current price is: \$3,151.00.



SKU: WINWHI

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Windsor toilet/commode lift seat provides a secure, safe, and user-friendly solution for transferring to the toilet, enabling individuals to use their home restrooms confidently. Utilizing the effective 'sit to stand' motion, this lift seat leverages the body's major muscle groups in the proper sequence, known as the most efficient way to rise to a standing position. It imitates the rise motion of standard "rise and recline" lift chairs, so it feels familiar to many users.

FEATURES

- Maintain a safe body position whilst using the toilet
- Cleaner and easier functionality, with remote control, included
- Comfortable, Stable lifting/lowering onto the toilet
- Reduced impact on joints
- Remote control for ease of use
- Removable portable bed pan for convenient bedside use
- Rear wheels for easy manoeuvrability (move from room to room if required)
- Increased personal confidence and prolongs your ability to live at home
- The battery-powered device that enables users with restricted mobility to get on and off the toilet at the touch of a button
- SOS Alarm Button

SPECIFICATIONS

Dimensions (L x W x H)
Weight

60.5 x 52.5 x 79.5cm

22kg



Max. Loading Weight
Distance between handrails
Support times for battery full charge
Working Voltage
Shape

200kg50.6cm150110-240V Adaptive Voltage

Support 33 degree lift and down

BROCHURE

Windsor - Toilet Commode Seat Information