



WINDSOR PORTABLE TOILET/COMMUNE LIFT SEAT

\$3,501.00



SKU: WINWHI

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Windsor toilet/commode lift seat provides a secure, safe, and user-friendly solution for transferring to the toilet, enabling individuals to use their home restrooms confidently. Utilizing the effective 'sit to stand' motion, this lift seat leverages the body's major muscle groups in the proper sequence, known as the most efficient way to rise to a standing position. It imitates the rise motion of standard "rise and recline" lift chairs, so it feels familiar to many users.

FEATURES

- Maintain a safe body position whilst using the toilet
- Cleaner and easier functionality, with remote control, included
- Comfortable, Stable lifting/lowering onto the toilet
- Reduced impact on joints
- Remote control for ease of use
- Removable portable bed pan for convenient bedside use
- Rear wheels for easy manoeuvrability (move from room to room if required)
- Increased personal confidence and prolongs your ability to live at home
- The battery-powered device that enables users with restricted mobility to get on and off the toilet at the touch of a button
- SOS Alarm Button

SPECIFICATIONS

Dimensions (L x W x H)

60.5 x 52.5 x 79.5cm

Weight

22kg

Max. Loading Weight	200kg
Distance between handrails	50.6cm
Support times for battery full charge	150
Working Voltage	110-240V Adaptive Voltage
Shape	Support 33 degree lift and down

BROCHURE

[Windsor - Toilet Commode Seat Information](#)