




WINNCARE PURE AIR 120 BARIATRIC MATTRESS



SKU: SM158ADU

PRODUCT DESCRIPTION

Purpose & Benefits

The **Pure Air 120 Bariatric Mattress** (200 cm × 135 cm) is expertly designed to deliver premium pressure relief and therapeutic support for plus-size patients at very high risk of pressure injuries. Featuring full-width alternating air cells with a pioneering figure-of-eight i-beam pattern, it provides smooth surface movement, optimal immersion, and even pressure redistribution. The whisper-quiet control unit, durable PU cell construction, and easy-clean, high-stretch cover ensure reliability and ease of use—making it ideal for high-risk patients and trusted by caregivers alike.

Features

- Full-width alternating air cells with figure-of-eight i-beam pattern to maximise pressure redistribution
- Therapy modes include alternating, static (auto-return), seat mode, and CPR rapid deflation
- Whisper-quiet pump with audible alerts: low pressure, cycle fault, power failure; plus control-panel lockout
- Welded, multi-stretch, waterproof and vapour-permeable cover; machine washable up to 95 °C
- External, easy-replace control-unit filters support in-house servicing

- Durable construction with a 3-year warranty (subject to annual servicing)

Specifications

- **Maximum User Weight:** Up to 298 kg
- **Dimensions (inflated):** 200 cm (L) × 135 cm (W) × 20 cm (H)
- **Alternate Width Options:** Also made to order in 120 cm (standard) or 150 cm widths
- **Warranty:** 3 years (with annual servicing)