

**WINNCARE PURE AIR 120 BARIATRIC MATTRESS -  
200CM (L) X 120CM****SKU:** SM158ADU**PRODUCT DESCRIPTION****Purpose & Benefits**

The **Pure Air 120 Bariatric Mattress** is expertly designed to deliver premium pressure relief and therapeutic support for plus-size patients at very high risk of pressure injuries. Featuring full-width alternating air cells with a pioneering figure-of-eight i-beam pattern, it provides smooth surface movement, optimal immersion, and even pressure redistribution. The whisper-quiet control unit, durable PU cell construction, and easy-clean, high-stretch cover ensure reliability and ease of use—making it ideal for high-risk patients and trusted by caregivers alike.

**Features**

- Full-width alternating air cells with figure-of-eight i-beam pattern to maximise pressure redistribution
- Therapy modes include alternating, static (auto-return), seat mode, and CPR rapid deflation
- Whisper-quiet pump with audible alerts: low pressure, cycle fault, power failure; plus control-panel lockout
- Welded, multi-stretch, waterproof and vapour-permeable cover; machine washable up to 95 °C
- External, easy-replace control-unit filters support in-house servicing

- Durable construction with a 3-year warranty (subject to annual servicing)

## Specifications

- **Maximum User Weight:** Up to 298 kg
- **Available Dimensions (inflated):**
  - SKU: **SMO23 - 105AU** | 200cm (L) x 105cm (W) x 20cm (H)
  - SKU: **SM158AU** 200cm (L) x 120cm (W) x 20cm (H)
  - SKU: **SM158DAU** 200cm (L) x 135cm (W) x 20cm (H)
- **Alternate Width Options:** Also made to order in 120 cm (standard) or 150 cm widths
- **Warranty:** 3 years (with annual servicing)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)