



# WINNCARE PURE AIR 120 BARIATRIC MATTRESS - 200CM (L) X 120CM



SKU: SM158ADU

#### **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The **Pure Air 120 Bariatric Mattress** is expertly designed to deliver premium pressure relief and therapeutic support for plus-size patients at very high risk of pressure injuries. Featuring full-width alternating air cells with a pioneering figure-of-eight i-beam pattern, it provides smooth surface movement, optimal immersion, and even pressure redistribution. The whisper-quiet control unit, durable PU cell construction, and easy-clean, high-stretch cover ensure reliability and ease of use—making it ideal for high-risk patients and trusted by caregivers alike.

#### **Features**

- Full-width alternating air cells with figure-of-eight i-beam pattern to maximise pressure redistribution
- Therapy modes include alternating, static (auto-return), seat mode, and CPR rapid deflation
- Whisper-quiet pump with audible alerts: low pressure, cycle fault, power failure; plus controlpanel lockout
- Welded, multi-stretch, waterproof and vapour-permeable cover; machine washable up to 95 °C
- External, easy-replace control-unit filters support in-house servicing



• Durable construction with a 3-year warranty (subject to annual servicing)

## **Specifications**

- Maximum User Weight: Up to 298 kg
- Available Dimensions (inflated):
  - SKU: SMO23 105AU | 200cm (L) x 105cm (W) x 20cm (H)
  - SKU: SM158AU 200cm (L) x 120cm (W) x 20cm (H)
  - SKU: SM158DAU 200cm (L) x 135cm (W) x 20cm (H)
- Alternate Width Options: Also made to order in 120 cm (standard) or 150 cm widths
- Warranty: 3 years (with annual servicing)