




## WONDERSHEET +PLUS

**\$111.00 - \$164.00** Price range: \$111.00 through \$164.00



**SKU:** WS1DB35-H

## PRODUCT DESCRIPTION

## PURPOSE & BENEFITS

The Wonder Sheet +PLUS is a unique bed sheet designed to help people with mobility issues get in and out of bed more easily. It has a satin panel in the middle that allows you to slide and pivot your hips over more easily without having to lift them. The sides of the sheet are made of a cotton/poly blend, which provides additional stability and prevents the sheet from slipping.

## FEATURES

- Reduces pain and discomfort
- Makes turning over in bed easier
- Assists with getting in and out of bed
- Requires less staff assistance
- Satin panel for effortless sliding and repositioning without lifting
- Cotton/poly blend sides for additional stability and slip prevention
- Machine washable and wrinkle-resistant

## SPECIFICATIONS

- **Materials:** Satin panel (90% polyester, 10% spandex), cotton/poly blend sides (50% cotton, 50% polyester)
- **Dimensions:** Available in long single, double, and queen sizes.
- **Weight:** 1.3kg