

OSSUR FORMFIT WRIST & FOREARM BRACE



SKU: OSS617083-H

PRODUCT DESCRIPTION

WRIST BRACE PURPOSE & BENEFITS

The Össur Formfit® Wrist Brace & Forearm Brace provides immobilisation and protection for patients suffering from tendonitis, select soft tissue injuries (including mild to moderate sprains and strains) as well as for post-cast healing. The lightweight, breathable material and contour fit make this comfortable product an excellent choice.

KEY FEATURES

- Breathable material
- Contour Fit
- Contact Closure Straps
- Customisable, mouldable Aluminium Stays

SIZING

Wrist Circumference:

Xsmall	13 - 16cm
Small	16 - 18cm
Medium	18 -20cm
Large	20 - 23cm
Xlarge	23 - 25cm

**measure wrist circumference at the ulnar styloid process*

What is a forearm brace for?

The brace takes stress off the tendon at the spot it meets the bone, which is often the most tender spot. This enables the painful point to cease which will decrease the pain throughout the performance and empower you to make the throw to first or swing through.

Purchase this Wrist & Forearm Brace online at Medimart, Also check out our [shoulder brace](#) range.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au