



## ZOOM KENSINGTON TOILET/COMMODE LIFT SEAT

**\$3,089.00**



**SKU: KENWHI**

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Zoom Mobility Kensington toilet/commode offers a secure and reliable lifting mechanism, which can be effortlessly controlled remotely. This user-friendly toilet transfer solution utilizes a familiar rising motion found in most lift chairs. As a result, individuals can use the toilet with ease and independence, requiring less assistance and reducing the likelihood of accidents such as falls or becoming trapped.

### FEATURES

- Ergonomically designed to relieve the knee force required to transfer to and from a toilet seat
- Removable portable bed pan for convenient bedside use
- The battery-powered device that enables users with restricted mobility to get on and off the toilet at the touch of a button
- Slowly lifts the body and drops down for added safety
- Armrest support enables greater stability and safety when transferring for toileting
- Rear wheels for easy maneuverability (move from room to room if required)
- Maximum Load weight 200kg

### SPECIFICATIONS

<b>Dimensions (L x W x H)</b>	60.5 x 52.5 x 79.5cm
<b>Weight</b>	20kg
<b>Max. Loading Weight</b>	200kg
<b>Distance between handrails</b>	50.6cm
<b>Support times for battery full charge</b>	150

**Working Voltage**  
**Shape**

110-240V Adaptive Voltage  
33-degree lift up and down

**BROCHURE**

[Kensington - Toilet Commode Seat Information](#)