




ZOOM KENSINGTON TOILET/COMMODE LIFT SEAT

\$3,089.00



SKU: KENWHI

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Zoom Mobility Kensington toilet/commode offers a secure and reliable lifting mechanism, which can be effortlessly controlled remotely. This user-friendly toilet transfer solution utilizes a familiar rising motion found in most lift chairs. As a result, individuals can use the toilet with ease and independence, requiring less assistance and reducing the likelihood of accidents such as falls or becoming trapped.

FEATURES

- Ergonomically designed to relieve the knee force required to transfer to and from a toilet seat
- Removable portable bed pan for convenient bedside use
- The battery-powered device that enables users with restricted mobility to get on and off the toilet at the touch of a button
- Slowly lifts the body and drops down for added safety
- Armrest support enables greater stability and safety when transferring for toileting
- Rear wheels for easy maneuverability (move from room to room if required)
- Maximum Load weight 200kg

SPECIFICATIONS

Dimensions (L x W x H)	60.5 x 52.5 x 79.5cm
Weight	20kg
Max. Loading Weight	200kg
Distance between handrails	50.6cm
Support times for battery full charge	150

Working Voltage
Shape

110-240V Adaptive Voltage
33-degree lift up and down

BROCHURE

[Kensington - Toilet Commode Seat Information](#)