

icare

IC100

The IC100 Bed is a non-function hospital grade bed, with a static mattress platform. This bed provides a strong base for a partner or caregiver to sleep alongside their patient.



BED CODES & SIZES:



LONG SINGLE

Code: IC100LS

Mattress Surface Size: 2030 x 900

Base Size: 2150 x 960



KING SINGLE

Code: IC100KS

Mattress Surface Size: 2030 x 1070

Base Size: 2150 x 1120

COLOUR CHOICES:



Onyx



Stone


KEY FEATURES:



200kg Weight
Rating



Fixed Base
Height 37cm



Available with
Headboards



Lockable Heavy
Duty Castors



Strong Breathable
Mattress Platform



End And Side
Mattress Retainer
Bars



10 Year Warranty

TECHNICAL SPECIFICATIONS:

- » Safe working load - The safe working load is 200kg includes any model of Icare mattresses.
- » Fabric - Hard wearing commercial grade fabric. Onyx range (BLACK) 82% Polyester, 18% Olefin. Stone range (GREY) 70% Polyester, 30% Acrylic.
- » Mattress Surface - Strong metal mesh under mattress allowing increased airflow.
- » Wheels - large lockable heavy duty swivel castors on each corner.
- » Base fixed height - 37cm from floor to top of base.

ACCESSORIES AVAILABLE:

- » High Side Rail
- » Low Side Rail
- » Bed Stick
- » Over Bed Pole
- » Full Length Side Rails
- » Headboards