

# **ICARE MATTRESSES**

**VISCO ELASTIC MATTRESSES:** These mattresses is heat and pressure sensitive which means that as it gets warmer it becomes softer. So when there is a pressure point developing which usually gets hot and enflamed, this area of the mattress will become softer which allows the blood to flow freely again.

Suitable for all Icare beds



#### **KEY FEATURES:**



Visco Elastic Memory Foam



200kg Weight Limit



2 Way Stretch Zip off Washable Cover



10 Year Warranty

## **MATTRESS FEELS:**



IC15 Firm



IC20 Medium



IC25 Soft

**M2 MEDICAL MATTRESS:** The M2 mattress features firm edge support and high density foams. This is a high risk mattress, and suitable for the use on all Icare beds.



## **KEY FEATURES:**



15cm Mattress Depth



200kg Weight Limit



Fully Waterproof Cover



5 Year Warranty



LATEX MATTRESS: The Latex mattresses are a low risk pressure mattress but provide a unique bounce back support which makes it easy for transfers.

Suitable on any base platform.



## **KEY FEATURES:**











18cm Mattress Depth 200kg Weight Limit

Natural Latex

Quilted Ultra 5 Year Warranty

MATTRESS OVERLAYS: The range of Icare Pressure care overlays are designed to provide added pressure relief and comfort on an existing mattress. These overlays are heat and pressure sensitive causing them to soften around pressure points and inflamed areas. These overlays come in all the mattress sizes up to queen.



### **KEY FEATURES:**



5cm Overlay Thickness



Visco Elastic Memory Foam



Suitable for any Mattress Type



2 Way Stretch zip off washable cover



10 Year Warranty