

# EN9000 Endless Bed

*A bed with Endless possibilities*




*...the true  
floorline bed  
with no towers...*

**TRUE FLOORLINE BED - NO TOWERS**

**ONE BED - ENDLESS POSSIBILITIES**


**AUSTRALIAN DESIGNED & MADE**



Unique dual Safe Working  
Load of 200kg or 300kg  
- simply swap handsets




Universal solution for  
Low, High & Bariatric care



Low height of 100 and  
high height of 725mm



Chair position  
at the push of a button



No towers  
Up to 400mm shorter  
than other floor level beds



10 year frame warranty &  
5 year Linak warranty

*We are proud to introduce the  
first ever universal bed that suits  
low care, high care and bariatric  
care. Plus, this innovative three-  
in-one bed has a low height of  
100mm and NO TOWERS!*

*From care of frail elderly or  
dementia patients, through to  
acute and bariatric settings,  
the patent pending Australian  
made EN9000 Endless Bed has  
you covered. One bed, endless  
possibilities!*

**ALRICK**  
HEALTHCARE  
1800 282 001

# EN9000 Endless Bed

*A bed with endless possibilities*

## ENDLESS BED POSITIONS



High position




Medium position



Low position



Backrest function



Backrest & Kneebreak




Backrest & Kneebreak



Trendelenburg



Reverse Trendelenburg



Chair position

## SPECIFICATIONS

Overall dimensions	2080mm x 905mm	Single
	2080mm x 1090mm	King single
Sleep surface dimensions	2000mm x 890mm	Single
	2000mm x 1070mm	King single
SWL	200kg/300kg	
NET Bed Weight	100kg	
Low Position	100mm	
High Position	725mm	
Underbed Clearance	Clearview	
Angles	Back	62°
	Knee	35°
	Trendelenburg	15°
	Reverse	15°

Warranty	Frame	10 years
	Electric	5 years
Optional Features	Underbed Light	
	Battery Backup	
	Dual Braking	
	End Board Options	
Product Codes	Single	EN9000
	King Single	EN9000KS
	Double	EN9000D
	Queen	EN9000Q
	King	EN9000K