

PICASSO RECLINER

**Four motor lift with
headrest & lumbar
adjustment functions**

The Picasso is the
combination of style and
support with comfort



General Conditions that the Power Headrest and Lumbar adjustment can assist with:

- Kyphosis
- Lordosis
- Arthritis

Customer Feedback:

I came over to your shop yesterday to try out the latest chair to hit Australia and wow what a chair! There is nothing in your shop the comes close and I've tried them all. Having a fused neck, arthritis and just had two lower back surgeries I am needing support while sitting. I could never get a chair to support my lower back and neck until now because of the electric adjustable lumbar support and electric adjustable head rest. Gone are the days of trying to use pillows to get the support which could never be achieved. Even the way the chair lifts you has improved as you don't feel like you are going to land on your face.

Therapist Feedback:

I feel this would have a good application for kyphosed clients and offer adjustability not in the typical lift chairs ranges on offer. The contoured seat base is also a good feature, to assist with lower limb and pelvic positioning, acting as a lateral thigh guide to a degree. Although often client with a higher or more complex seating needs may not be the best candidate for the standard lift chairs ranges, these will be welcome additions to provide low to moderate seating support and correction which I think will be beneficial.

4 INDEPENDENT MOTORS FOR A WIDE RANGE OF POSITION COMBINATIONS AND SUPPORT

Power Headrest

The individually adjustable head support allows you to lift or lower the position of the headrest even when the backrest is fully reclined providing a higher level of adjustment for your neck and back.

Power Lumbar

Individually adjustable lumbar support, which can assist in maintaining a proper sitting posture; helping to relieve tension on the spine to prevent pain or discomfort due to poor posture, arthritis or existing back problems.

Power Backrest and Power Footrest

With a the dual lift system, it allows you to individually adjust the backrest and the footrest for the position that best fits the clinical needs of the client



The new four motor lift mechanism has a more vertical lift with a slight tilt compared to the forward tilt you find on generic chairs



FEATURES

- Dual lift function
- Lay-flat recliner
- Footrest extension
- All close & all lift buttons
- Hand controller
- Ultra durable fabrics
- Lumbar support to accommodate sitting posture

Easy to use hand control

- Soft grip buttons
- 6 pre-defined positions
- Choose “full recline” or “full stand” with the touch of a button.



SPECIFICATIONS

Description	Overall Depth	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Reclined Length	Product Weight	Max User Weight
Picasso Lift Recliner KA551	1003mm	813mm	1080mm	508mm	521mm	495mm	1830mm	74kg	158kg