

ICARE IC280

The IC280 Bed combines the essential features of a homecare bed, with ultra high lift function enabling easier patient care. The three function design is compatible with all Icare accessories and attachments.

KEY FUNCTIONS



Head Lift



Knee Break



Height Adjustable
38-81cm



Note: Standard Castors are 14cm,
or high option 16cm.

BED CODES & SIZES



LONG SINGLE
Code: IC280LS
Mattress Surface Size: 2030 x 900
Base Size: 2150 x 960



KING SINGLE
Code: IC280KS
Mattress Surface Size: 2030 x 1070
Base Size: 2150 x 1120


COLOUR CHOICES



ONYX



STONE



SMOKE



SILVER



VINYL EBONY



VINYL DOVE

KEY FEATURES



180kg Safe Working Load



Strong Breathable Mattress Platform



Available with Headboards & Footboards



External Rechargeable Battery Back-up



Heavy Duty Castors



Vertical Lift in Place



End and Corner Mattress Retainer Bars



5 Year Warranty



Back Lit Wipeable Hand Piece



No gap between partner mattresses

TECHNICAL SPECIFICATIONS

- Safe working load - 180kg on long single and king single sizes. The safe working load includes any model of Icare mattresses.
- Hand Piece - Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging to side of bed. Wipeable.
- Fabric - Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties.
- Mattress Surface - Strong metal mesh under mattress allowing increased airflow
- No pull motors for added safety.
- Wheels - 4x heavy duty swivel castors with two lockable castors.
- Head lift angle 62° and knee break angle 42°
- Max height - 78cm from floor to top of mattress platform.
- Min height - 35cm from floor to top of mattress platform.

ACCESSORIES AVAILABLE

- High Side Rail
- Low Side Rail
- Bed Stick
- U-Assist Rail
- Over Bed Pole
- Over Bed Table
- Full Length Side Rails
- Headboard
- Footboard