


## VERSAFLO3® SYSTEM




## Pressure Redistribution Information

## **VERSAFLO3® SYSTEM**

Risk Level: High SWL Minimum 20kg SWL Maximum 140kg

Waterlow ≤ 20 Immersion Rating: Medium Envelopment Score: 4/5



Risk category ratings, immersion ratings and envelopment scores area based on average peak pressure, therapist recommendations case studies and Waterlow scale. For some high risk clients Icare Zerotec covers are recommended. The best way to assess a client on a pressure reducing surface is to monitor pressure points with visual checks over a period of five days.