Royale Flush

LUXURY TOILET BIDET LIFT SEAT



Top Gun Mobility
Royale Flush
toilet/commode lift
seat delivers a stable,
safe, easy to use toilet
transfer solution that
restores an individual's ability to
safely toilet in the comfort of the

safely toilet in the comfort of their home.
This Electronic Toilet seat riser replicates natural standing movement at the touch of a button, assisting users to help themselves get up into a standing position.
The Royale Flush comes with built in Bidet for that extra level of comfort and luxury.



Royale Flush Features

- Safe toileting at home giving you independence
- Lowers you under power from a standing position until you are positioned over the toilet
- Raises you back up under power from a sitting position
- Safely use the toilet with reduced assistance or on your own
- Helps reduce manual handling risks for carers who provide toileting assistance
- Reduce the risk of falls and associated injuries, and reduce the likelihood of be unable to get up from the toilet
- Electric washing Bidet feature included (plumbing required)
- Remote control for ease of use
- Maximum Load weight 200kg

Specifications

Dimensions (L x W x H)	60.7 x 61.6 x 70.6cm
Weight	22kg
Max. Loading Weight	200kg
Distance between handrails	50.4cm
Support times for battery full charge	N/A Unit is plugged into power outlet
Working Voltage	240V adaptive voltage
Shape	33 degree up and down

NOTE: All specifications subject to change without notice.

- The information contained herein is correct at the time of publication; due to our commitment to constantly improve our products during development, we reserve the right to alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.

