



# Premium rollators



EN

made in Norway



# We are proud to present our revitalized and modern rollator portfolio

The updated versions the TOPRO Troja 5G and TOPRO Olympos ATR, which replaces the TOPRO Troja 2G and TOPRO Olympos, demonstrate that TOPRO is consolidating its position as market leader in the premium rollator segment.

We are continuing our product updating programme with the new, modern designs which we began with TOPRO Pegasus in the autumn of 2019. We are enthusiastically introducing many updated practical features that our customers have requested.

TOPRO has, for more than 30 years, worked with the development, production, and sale of quality mobility aids. We are proud of the products we make, and we are happy to see the benefits these products provide for our customers.

Through a continuous focus on improvements for the benefit of our customers and the environment, we have also placed significant emphasis on sustainability in our production process. Our factory in Norway uses hydropower and environmentally friendly aluminium to minimize the carbon footprint. Our injection moulding plant uses clean materials that are recyclable to contribute

to the circular economy. We have also implemented a “zero waste” principle to ensure no plastic contaminates the natural environment for the benefit of future generations.

We hope that our products contribute to stable, safe walking and that they motivate everyone to „Keep on moving“.



**Asbjørn Eskild**  
CEO TOPRO Industri AS



# Advantages of the TOPRO rollator family

Our mobility and design specialists are constantly working to improve our products with optimized functionality and modern designs for the twenty twenties.

We believe that good design, new colour options, practical high-quality functions, give added value to our products.

Safety and functionality are our top priorities and we want you to feel safe and comfortable using a TOPRO rollator for years to come.

## Unique memory function

Saves the preferred preset handle height, for easy re-use. (Standard on Troja 5G, Olympos ATR and Troja Classic)

## Easy wheel change

The interior braking system (TOPRO IBS) in combination with the TOPRO Quick Release system for easy wheel change, allows the use of different tyre treads. (Standard on TOPRO Pegasus, Troja 5G, Olympos ATR and Odysse)

## Ergonomic handles

The TOPRO Ergo Grip ensures correct walking posture and comfortable armrests while seated. Together with the soft tyres it provides double shock absorbing effect for stable driving comfort. (Standard on TOPRO Pegasus, Troja 5G, Olympos ATR and Odysse)

## Easy to fold

Simply pull the strap on the seat with only one hand. Practical to transport, stays firmly on its wheels and takes up little space when stowed.

## Tilt function and edge guard

Enables the user to easily negotiate thresholds and kerbs, while keeping both hands on the rollator for safety.



# TOPRO PEGASUS

CARBON FOR LIGHTWEIGHT, TOPRO FOR QUALITY



Modern, slender and lightweight carbon rollator with functionality and elegant design. This design combines to guarantee the highest level of quality, safety and stability at a weight of just 6.1 kg\*.

Tried-and-tested features have been further developed:

- User friendly method for setting the height of the handles using a **numeric scale**.
- **Partly integrated brake cables** and modern reflective panels on the frame.
- Elegant detachable shopping bag, which, with its retracted position provides **increased forward view** to **improve user safety**.
- Easy handling and manoeuvrability.
- **Reflectors** on the frame, lock wheels and shopping bag provide **360 degree visibility for safety in the dark**.
- **Shock-absorbing tyres** and **ball bearings** in wheels and forks for a comfortable ride.

Available in colour:



White



Steel blue

Measurements, colour options and other product details are available on page 19.





“TOPRO Pegasus represents a new trend in the development of mobility aids where design, functionality and safety are combined to meet users’ expectations.”



made in Norway





### **Elegant and practical shopping bag**

The elegant rollator design continues with the new shopping bag. Smart and practical, with internal pocket providing a storage option for small items, such as keys and a purse. Carries up to 10 kg.



# New features of the TOPRO Pegasus

## Safety and quality

Safety, high durability, a patent-pending brake and various sophisticated functions are notable features of the TOPRO Pegasus. The modern, slender and lightweight design makes it unique.



Numeric scale enables user friendly height adjustment of the handles.





# TROJA 5G

THE COLOURFUL ALLROUNDER



Upgraded level of handling comfort and stability of the ground breaking Troja Classic. The TOPRO Troja 5G offers a modern, slender design and a high level of quality, safety and driving comfort.

Innovative upgrades to provide increased comfort and ease of use:

- User friendly method for setting the height of the handles using a **numeric scale**.
- **Shock-absorbing tyres** provides a comfortable ride and smooth suspension.
- **Bell** to attract attention when necessary. Easy and safe to operate from the handle.
- Easy handling and manoeuvrability due to **ball bearings in forks and wheels**.
- **Detachable shopping bag** with reflective band, zipper and shoulder strap. Closeable and holds up to 10 kg. Two inner pockets for small items.
- Easy to fold, stays firmly on its wheels and takes up little space when stored.

Measurements, colour options and other product details are available on page 20.

Model M & S available in on-trend matt coated colours:





“I have had a TOPRO rollator for a year now and it has transformed my life. I used to have a stick before, but once I got the rollator, I found I could walk so much better.”

| Audrey from the UK



made in Norway



# The new TOPRO Troja 5G



The slender design gives a modern and light appearance – matching it perfectly to the revitalized TOPRO rollator family.

We know how important it is to adapt the rollator to individual style and needs. With the new product line, we are proud to offer a wide range of colour options and accessories to help you find your favourite TOPRO rollator and keep on moving.



The detachable shopping bag is redesigned with an on-trend grey fabric print on the lid. Reflector band for visible safety in the dark and practical carrying strap.

New on-trend matt colour coating and a wide range of colour\* options available:





Numeric scale for easy height adjustment of the handles, with improved lock clamp.

To ensure that our users are comfortable walking with a rollator, TOPRO's rollator specialists have given it an attractive design. The new frames have a more slender and lighter appearance and the new wheels a slimmer design, while retaining TOPRO's well-known stability and shock-absorbing function.





# TROJA CLASSIC

THE ORIGINAL



The Original and first innovative rollator on the market to offer enhanced functionality and design that folds sideways. Is still as popular as ever and a top-selling premium rollator for over 15 years.

The design from TOPRO that inspired a whole new generation of rollators:

- **Lightweight** at only 6.8 kg, and with a **weight bearing capacity** of 165 kg.
- Innovative **sideways folding technology** that is easy to operate with one hand by pulling the strap on the seat.
- **Ball bearings in forks and wheels** to ensure an easy ride and manoeuvrability.
- **Practical net basket** for easy access to items with **reflective band**.
- Available in three different sizes.
- Valued by users for being easy to handle and operate, having great manoeuvrability which makes it a trusty companion in everyday situations.
- It is also valued for consistently performing well in tests and being a previous test winner\*.

**Model M available in:**



**Model S, X available in:**



Measurements, colour options and other product details are available on page 21.



\* Test 1; Conducted by the German consumer organisation Stiftung Warentest (9/2005).  
Test 2; Conducted by the German consumer programme Marktcheck in 2014, having previously been recognised by the magazine Ökotest and the specialist gerontotechnology company GGT®.



made in Norway



# OLYMPOS ATR

THE ALL-TERRAIN ROLLATOR



With large front wheels the Olympos ATR provides reliable and comfortable support on uneven roads and cobblestones as well as on sand or gravel foundations. The perfect choice for outdoor activities.

Innovative features to provide maximum support for outdoor use:

- **Wider design** to ensure stability in an outdoor environment.
- **Large and shock-absorbing tyres** to enable a comfortable ride and maintain stable and smooth manoeuvrability on rough terrain.
- **Detachable spacious water-resistant shopping bag** in on-trend grey fabric print. With reflective band and zipper. Closeable and holds up to 15 kg.
- User friendly method for setting the height of the handles with **numeric scale**.
- Easy to ride and manoeuvre due to **ball bearings in forks and wheels**.

Measurements, colour options and other product details are available on page 22.

Model M & S available in on-trend matt coated colours:



Brand new on the Olympos ATR – the ergonomic handle (**TOPRO Ergo Grip**) makes it easier and safer to get up from the sitting position. Together with the soft tyres it provides double shock absorbing effect for stable driving comfort.





made in Norway



# ODYSSÉ

THE TRAVEL ROLLATOR



The Odyssé be takes up no more space than a sports bag when folded for transport, can easily fit into a car boot or stowed on a coach or train and easily reassembled in seconds!

Innovative features which make the TOPRO Odyssé the ideal travel companion:

- **Compact size** for transport with its patented 2-way folding system that folds the rollator both vertically and horizontally.
- Practical **transport bag** included as standard.
- Can be quickly unfolded when needed, and is **small and easy to handle** for storage.
- Fits in the ordinary luggage compartment on trains so does not block the aisles.
- Due to its compact size some airlines allow it to be stored as hand luggage in the cabin. Remember to ask the airline in advance!
- **Spacious, expandable shopping bag** in water-resistant fabric. Large cover with **reflective band**. Holds up to 10 kg.
- Easy to push and manoeuvre due to **ball bearings in forks and wheels**.

---

Measurements, colour options and other product details are available on page 23.

Available in colour:



Ultramarine



Wine red





“The TOPRO Odyssé  
is a great travel companion.”



made in Norway





# Accessories



**Crutch holder**  
can be fitted on either side.

TOPRO Pegasus Art. No. 815236  
Troja Classic Art. No. 814025  
Troja 5G/  
Olympos ATR Art. No. 814040  
Odyssé Art. No. 814457

TOPRO Pegasus, Troja 5G,  
Olympos ATR, Troja Classic,  
Odyssé, Troja Neuro



**Tray**  
Snaps in securely, easy to remove  
and clean, carries up to 5 kg

TOPRO Pegasus, Troja Classic/5G,  
Odyssé Art. No. 814728  
Olympos ATR M Art. No. 815212

TOPRO Pegasus, Troja 5G,  
Olympos ATR, Troja Classic,  
Odyssé



**Back support with padding**  
ergonomically height adjustable  
according to preferred seating  
position. Available in two  
lengths.

short 67 cm Art. No. 815270  
long 75 cm Art. No. 815272

TOPRO Pegasus, Troja 5G, Olympos  
ATR, Odyssé, Troja Neuro



**LED lamp with clip on and USB  
plug for safety in the dark.**  
Chargeable.

Art. No. 815260

TOPRO Pegasus, Troja 5G,  
Olympos ATR, Troja Classic,  
Odyssé, Troja Neuro



**Back support with padding**  
ergonomically height adjustable  
according to preferred seating  
position. Available in two  
lengths.

short 68 cm Art. No. 815271  
long 76 cm Art. No. 815273

Troja Classic



**Umbrella with attachment**  
keep dry under the umbrella  
in rainy weather or use it as a  
parasol.

Art. No. 814068

TOPRO Pegasus, Troja 5G,  
Olympos ATR, Troja Classic,  
Odyssé, Troja Neuro

Additional accessories can be found on our website.

# Technical data for TOPRO Pegasus



<b>Recommended user height</b>	150–200 cm
<b>Max. user weight</b>	130 kg
<b>Max. capacity shopping bag</b>	10 kg / 15 l
<b>Handle height</b>	79–97.6 cm
<b>Width between handles</b>	47 cm
<b>Handles</b>	Ergonomic, several holding positions (TOPRO Ergo Grip)
<b>Maximum length</b>	65 cm
<b>Maximum width</b>	64 cm
<b>Seating width</b>	43.5 cm
<b>Seat height</b>	63 cm
<b>Turning circle</b>	84 cm
<b>Diameter wheels</b>	20 cm
<b>Width wheels</b>	3.6 cm
<b>Grip distance driving brake</b>	7 cm
<b>Height when folded</b>	80.3 cm
<b>Width when folded</b>	27.3 cm
<b>Length when folded</b>	65 cm
<b>Transport dimensions (l × w × h)</b>	65 × 27.3 × 80.3 cm
<b>Product weight (excl. bag)</b>	6.1 kg
<b>Product weight (incl. bag)</b>	6.6 kg
<b>Area of use</b>	Indoors & outdoors
<b>Article No.</b>	815020 + colour code

**Available in colour:**



White



Steel blue

# Technical data for Troja 5G



Size	M (medium)	S (small)
Recommended user height	150–200 cm	135–170 cm
Max. user weight	150 kg	125 kg
Max. capacity shopping bag	10 kg / 15 l	10 kg / 15 l
Handle height	78–100 cm	71–86 cm
Width between handles	40 cm	40 cm
Handles	Ergonomic, serveral holding positions (TOPRO Ergo Grip)	
Maximum length	65 cm	65 cm
Maximum width	60 cm	60 cm
Seating width	45.5 cm	45.5 cm
Seat height	63,5 cm	54 cm
Turning circle	84 cm	84 cm
Diameter wheels	20 cm	20 cm
Width wheels	3.6 cm	3.6 cm
Grip distance driving brake	7 cm	7 cm
Height when folded	81.5 cm	74 cm
Width when folded	24 cm	24 cm
Length when folded	65 cm	65 cm
Transport dimensions (l × w × h) (excl. wheels and shopping bag)	60 × 24 × 76 cm	60 × 24 × 72 cm
Transport weight (excl. wheels and shopping bag)	5.6 kg	5.4 kg
Product weight (incl. wheels and shopping bag)	7.8 kg	7.6 kg
Area of use	Indoors & outdoors	
Article No.	814600 + colour code	814601 + colour code

Model M & S available  
in on-trend matt coated  
colours:





# Technical data for Troja Classic



Size	M (medium)	S (small)	X (extra small)
<b>Recommended user height</b>	150–200 cm	135–170 cm	125–160 cm
<b>Max. user weight</b>	165 kg	125 kg	80 kg
<b>Max. capacity basket</b>	10 kg / 20 l	10 kg / 20 l	5 kg / 11 l
<b>Handle height</b>	78–100 cm	71–86 cm	67–79 cm
<b>Width between handles</b>	47 cm	47 cm	37 cm
<b>Handles</b>	Anatomical		
<b>Maximum length</b>	65 cm	65 cm	65 cm
<b>Maximum width</b>	60 cm	60 cm	52 cm
<b>Seating width</b>	45.5 cm	45.5 cm	36 cm
<b>Seat height</b>	62 cm	54 cm	50 cm
<b>Turning circle</b>	84 cm	84 cm	79 cm
<b>Diameter wheels</b>	20 cm	20 cm	20 cm
<b>Width wheels</b>	3.6 cm	3.6 cm	3.6 cm
<b>Grip distance driving brake</b>	6.5 cm	6.5 cm	6.5 cm
<b>Height when folded</b>	80 cm	74 cm	65 cm
<b>Width when folded</b>	24 cm	24 cm	24 cm
<b>Length when folded</b>	65 cm	65 cm	65 cm
<b>Transport dimensions (l × w × h)</b>	65 × 24 × 80 cm	65 × 24 × 74 cm	65 × 24 × 65 cm
<b>Product weight (incl. wheels and basket)</b>	6.8 kg	6.6 kg	6.1 kg
<b>Area of use</b>	Indoors & outdoors		
<b>Article No.</b>	814750 + colour code	814770 + colour code	814740

## Model M available in:



Dark grey



Wine red



Silver

## Model S, X available in:



Silver

# Technical data for Olympos ATR



Size	M (medium)	S (small)
<b>Recommended user height</b>	150–195 cm	140–170 cm
<b>Max. user weight</b>	150 kg	150 kg
<b>Max. capacity shopping bag</b>	15 kg / 12 l	15 kg / 10 l
<b>Handle height</b>	78.5–94.5 cm	70.5–86.5 cm
<b>Width between handles</b>	42 cm	43 cm
<b>Handles</b>	Ergonomic, several holding positions (TOPRO Ergo Grip)	
<b>Maximum length</b>	77 cm	77 cm
<b>Maximum width</b>	68 cm	68 cm
<b>Seating width</b>	47 cm	47 cm
<b>Seat height</b>	61 cm	54.5 cm
<b>Turning circle</b>	97 cm	97 cm
<b>Diameter wheels</b>	28 cm / 23 cm	28 cm / 23 cm
<b>Width wheels</b>	4.2 cm	4.2 cm
<b>Grip distance driving brake</b>	7 cm	7 cm
<b>Height when folded</b>	87.5 cm	85 cm
<b>Width when folded</b>	27 cm	27 cm
<b>Length when folded</b>	77 cm	77 cm
<b>Transport dimensions (l × w × h)</b>	77 × 27 × 87.5 cm	77 × 27 × 85 cm
<b>Product weight</b> (incl. wheels and shopping bag)	8.9 kg	8.7 kg
<b>Area of use</b>	Outdoors	
<b>Article No.</b>	814300 + colour code	814305 + colour code

**Model M & S available  
in on-trend matt coated  
colours:**



# Technical data for Odyssé



<b>Recommended user height</b>	150–200 cm
<b>Max. user weight</b>	150 kg
<b>Volume capacity of shopping bag</b>	13–22.5 l
<b>Max. load in shopping bag</b>	10 kg
<b>Handle height</b>	81.5–100 cm
<b>Width between handles</b>	47 cm
<b>Handles</b>	Ergonomic, several holding positions (TOPRO Ergo Grip)
<b>Maximum length</b>	65 cm
<b>Maximum width</b>	60 cm
<b>Seating width</b>	45.5 cm
<b>Seat height</b>	62 cm
<b>Turning circle</b>	84 cm
<b>Diameter wheels</b>	20 cm
<b>Width wheels</b>	3.6 cm
<b>Grip distance driving brake</b>	7 cm
<b>Height when laterally folded</b>	81.5 cm
<b>Transport dimensions completely folded (l × w × h)</b>	75 × 24.5 × 44.5 cm
<b>Transport weight</b> (excl. wheels and shopping bag)	6.0 kg
<b>Product weight</b> (incl. wheels and shopping bag)	8.2 kg
<b>Area of use</b>	Indoors & outdoors
<b>Article No. Ultramarine / Wine red</b>	815000 / 815002

**Available in colour:**



Ultramarine



Wine red



## Premium rollators from Norway

### Manufacturer:



TOPRO Industri AS

Rambekkevegen 5

2816 Gjøvik

Norway

Tel. (+47) 61 13 46 00

info@topromobility.com

www.topromobility.com



Norwegian design and made in Norway.

TOPRO vouches for top quality and offers a 7-year warranty (excluding assembly).

Setup and variations in fixed accessories for the rollators may differ between countries.

### Quality Assurance Standards

EN ISO 9001 – quality

EN ISO 14001 – environment

EN ISO 11199-2 – rollator

EN ISO 10993-5 – biocompatibility

made in Norway



Dealer