



# Ultra Blood Pressure Monitor

ES-625B

Smart blood pressure measurement  
and pulse rate measurement.

- 5 Year Warranty\*
- Guest User Option
- Irregular Heartbeat Detection
- Automated Cuff Wrap Detection
- Measuring While Inflating



## SPECIFICATIONS:

**Power supply:** Battery: 6V DC 4x AAA batteries, AC adapter powered mode: 5V=1A (not included)

**Display mode:** Digital LCD Display V.A. 46 mm x 57 mm

**Measurement mode:** Oscillographic testing mode

**Measurement range:** Rated cuff pressure: 0 mmHg ~ 299 mmHg

**Measurement pressure:** SYS: 60 mmHg ~ 230 mmHg  
DIA: 40 mmHg ~ 130 mmHg

**Pulse value:** (40-199) beat/minute

**Accuracy:** Pressure: 5°C ~ 40°C within  $\pm 3$  mmHg, pulse value:  $\pm 5\%$

**Operating Condition:** A temperature range of: +5°C to +40°C  
A relative humidity range of 15% to 90%, non-condensing, but not requiring a water vapour partial pressure greater than 50 hPa

**An atmospheric pressure range of:** 700 hPa to 1060 hPa

**Storage & Transportation Condition:** Temperature: -20°C to +60°C  
A relative humidity range of  $\leq 93\%$ , non-condensing, at a water vapour pressure up to 50 hPa

**Measurement perimeter of the upper arm:** About 22-32 cm or 22-42 cm



Battery  
indication



Auto power  
on and off



Irregular  
heartbeat alert



Pulse rate  
display



250 memory  
storage

## AUSTRALIA

email: support@mylifsmart.net.au

Ph: 1300 110 155

(AEST - Mon-Fri, 11 am - 4 pm)

www.mylifsmart.net.au

## NEW ZEALAND

email: support@mylifsmart.net.nz

Ph: 0800 100 999

(NZST - Mon-Fri, 11 am - 4 pm)

www.mylifsmart.net.nz

## INDIA

email: support@mylifsmart.net

www.mylifsmart.net