



P001G

## Power Pack for Manual Chair

 **merits**  
**AUSTRALIA**

[www.meritsaustralia.com](http://www.meritsaustralia.com)



P001G

# Power Pack for Manual Chair

LED battery level indicator, speed adjustable knob, key switch, and hand throttle lever.



## Features

- Easy to attach and remove to most manual wheelchairs to allow attendant power assistance.
- Convenient forward & reverse switch & speed control, with easy to operate thumb throttle.
- Portable & lightweight for easy transportation which enables extra mobility anywhere.
- Wider twin wheels which allows much greater traction on grass and gravel.
- Universal mounting kit available as an option

## Technical Specs

Max Speed	5km/h / 3mph
Range	16km / 10mi.
Weight Capacity	165kg / 364lb
Item Weight	17kg / 37lbs
Motor	DC 12V 120W
Battery	12V / 22AH
Charger	3A off-board
Drive Wheel	8"x3" PU tire x 2 pcs