

Alrick EN9 Floorline Bed

Product Code: 7-EN9

The Alrick EN9 Floorline Bed is the first ever universal bed that suits low care, high care and bariatric care. This innovative three-in one bed has a low height of 100mm! One bed, endless possibilities.



Scan to see
more details



Product Features:

- Australian Made and Owned
- Floorline Model
- Vertical Lift
- Electric Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- SWL 200 - 250kg
- Back Lit Handset
- Slender Arch Headboard & Square Footboard
- Ironstone Powder Coated Frame



Forward & Reverse
Trendelenburg



Chair
Position



75mm Twin
Wheel Braking
Castors



Bluetooth
Connectivity
***handset must be
ordered separately***


Optional Features:

- CPR Backrest
- Electric Vascular Support
- Dual Locking Castors
- Bluetooth Handset (HC40)
- Attendant Control
- Battery Back Up
- Under Bed Light (Warm White)
- Underbed Light with Custom Colour
- USB Charging Ports
- Custom Powdercoat Colour
- Custom Endboard Style and Colour
- Luxa Full Surround (Timber or Upholstered)
- Floorline Partner Bed Set Up
- Custom Length Sleeping Surface (2m+)




Stocked Endboard Styles & Colours

Custom Endboard styles & colours available upon request.



Laminate
Slender Arch Style



Upholstered
Jersey Headboard &
Plain Footboard




Timber
Eltham Style



Natural Oak
Laminate



Ocean
Upholstered




Ash
Upholstered



Burnt Walnut
Timber



Beech
Timber



Saddle
Timber



Accessories – Added after manufacture

- 300kg SWL Handset (Single and King Single Only)
- Bluetooth Handset (7-BT-4-HB200)
- Mattresses
- Horizontal 400HR EN Rails (pair)
- Padded Cotside Protectors (pair)
- EN9 StandEzy Hook Handle
- Self Help Pole
- 200mm Bed Extension
- 200mm Foam Bolster
- Emergency Power Bed Unit
- IV Pole
- Combo Grab Handle
- EN T-Bar Buffer
- Reading Lamp



Scan to see
more details

Specifications:

Size	Single	King Single	Extra Wide	Double	Queen Double	King Double
External Dimensions	920 x 2115 mm	1090 x 2115 mm	1220 x 2115 mm	1390 x 2115 mm	1390 x 2115 mm	1820 x 2115 mm
Sleeping Surface	900 x 2000 mm	1070 x 2000 mm	1200 x 2000 mm	1370 x 2000 mm	1370 x 2000 mm	1800 x 2000 mm
Hi-Lo Position	100 – 725 mm	100 – 725 mm	170 – 725 mm	170 – 725 mm	170 – 725 mm	170 – 725 mm
Backrest / Knee Angle	62° / 35°	62° / 35°	62° / 35°	62° / 35°	62° / 35°	62° / 35°
TR / ATR Angle	15° / 15°	15° / 15°	15° / 15°	15° / 15°	15° / 15°	15° / 15°
SWL	200 kg	200 kg	250 kg	250 kg	250 kg	250 kg