



medimart
medical | mobility | therapy

Caremed Alrick.

Medimart Brighton

343 Brighton Rd, North Brighton
SA 5048

T: 08 7079 8773

E: brighton@medimart.com.au

Medimart Elizabeth Vale

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale, SA 5112

T: 08 7079 8777

E: elizabethvale@medimart.com

Medimart Reynella

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

T: 08 7079 8771

E: reynella@medimart.com.au

Medimart Monash

23 Winterton Road, Clayton, VIC 3168

T: 03 7503 8354

E: monash@medimart.com.au

Medimart Frankston

111 Cranbourne Rd, Frankston VIC 3199

T: 03 8763 3451

E: frankston@medimart.com.au

For further information go to
medimart.com.au
or Give us a call: **1300 633 063**

Alrick EN7
Endless Bed Series

‘True floorline bed
with no towers and
endless possibilities’



The EN7 series bed is a variant of the EN9000 manufactured overseas. It is the first ever universal bed that suits low care, high care and bariatric care, plus the innovative three-in-one bed. The EN7 uses Linak motors and has a low height of 120mm and has no towers!

Configuration Possibilities



Features

Electric Hi/Low	Auto regression backrest
Electric kneebend	75mm Twin wheel castors
Electric forward & reverse Trendelenburg	Standard head and foot boards (Fini Blonde laminate)

Specifications

Overall dimensions (mm)	2080 x 920 Single	
	2080 x 1090 King single	
Sleeping surface dimensions (mm)	2000 x 890 Single	
	2000 x 1070 King single	
SWL	200kg/300kg	
NET Bed weight	110kg	
Low Position	120mm	
High Position	725mm	
Underbed clearance	Clearview	
	Angles (°)	Back 62
		Knee 35
		Trendelenberg 15
		Reverse 15
Warranty	Frame	10 yr
	Electrics	5 yr
	Handset & Castors	1 yr
Optional features	Underbed light	
	Battery backup	
	Dual braking	
	End board options	

- 01 Electric forward & reverse Trendelenburg
- 02 Electric Kneebend
- 03 Auto regression backrest
- 04 Standard head & foot boards
- 05 10 year frame & 5 year electrical warranty
- 06 Electric Hi-Low
- 07 75mm Twin wheel castors

