

Presenting the Cobalt Shower Commode

COBALT

 HEALTH

 HELPING
 CARE
 HAPPEN

Attendant propelled, stainless steel shower commode

The Cobalt mobile shower commode is designed to provide a safer, more independent and more comfortable way to shower, transfer, and use the toilet.

Features

Versatile enough to be used in the shower, over the toilet or as a bedside commode

Easy to manoeuvre between bedroom and bathroom

Choice of 2 seat widths 46cm or 55cm

Anti-bacterial seat, backrest and handgrips

Soft, removable, polyurethane open front seat

Slide out, weight bearing foot plate


Quality stainless steel construction

4 x brake castors for added stability

Armrests flip up, and handgrips can be adjusted upwards, downwards or inwards, for comfort and safety

Warranty: Frame 5 years, castors, seat, handgrips, and back rest straps 1 year

SWL 200kgs – 46cm Seat, SWL 300kgs – 55cm Seat



Pan and
Pan Carrier
Included

Codes

46cm SEAT CODE: D01-SS-COM-46

55cm SEAT CODE: D01-SS-COM-55

Specifications

Code	Seat Width	Between Arms	Seat Height	Overall Width	Overall Height	Backrest Height	Unit Weight	SWL
D01-SS-COM-46	46cm	49cm	58cm	53cm	102cm	44cm	15kg	200kg
D01-SS-COM-55	55cm	58cm	58cm	62cm	102cm	44cm	17kg	300kg

Medimart Brighton

343 Brighton Rd, North Brighton
SA 5048

T: 08 7079 8773

E: brighton@medimart.com.au

Medimart Monash

23 Winterton Road, Clayton, VIC 3168

T: 03 7503 8354

E: monash@medimart.com.au

Medimart Elizabeth Vale

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale, SA 5112

T: 08 7079 8777

E: elizabethvale@medimart.com

Medimart Frankston

111 Cranbourne Rd, Frankston VIC 3199

T: 03 8763 3451

E: frankston@medimart.com.au

Medimart Reynella

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

T: 08 7079 8771

E: reynella@medimart.com.au

For further information go to

medimart.com.au

or Give us a call: **1300 633 063**