

AIRCOMFORT DELUXE BED STANDARD - ELECTRIC - INCLUDING POWER DRIVE

SKU: AC59112PD

PRODUCT DESCRIPTION

Purpose & Benefits

The **Aircomfort Deluxe Bed Standard - Electric - Including Power Drive** is designed to provide advanced pressure care and comfort for patients. Using its patented Airgonomic Seating System, it immerses and cushions the user, while powered adjustments and specialised materials support safe positioning, easier transfers and improved skin protection.

Features

- Patented Airgonomic Seating System with four air cells for immersion and cushioning
- Electric adjustments for back recline, seat tilt and raising legrest
- Drop-away and adjustable armrests and back wings for easier transfers
- VersoTela Fusion material designed for pressure care and sensitive skin
- Large castors for easier manoeuvrability
- One directional lock castor assists with straight line tracking
- Customisation and colour options available

Specifications

- **Drive Direction:** Forwards & backwards
- **Max Assist Speed:** Up to 6 km/h
- **Max Advised Ramp Angle:** 6°Wheel Lift Speed: 4 seconds up/down
- **Product Weight (excl. charger):** 18 kg
- **Force of Wheel on Floor:** 40 kg
- **Battery Voltage / Capacity:** 24V DC / 9Ah NiMH
- **Battery Lifetime:** 500 cycles or ~2 years
- **Average Running Distance:** ~9 km on full charge (varies by load, slope & use)
- **Safety Features:** Emergency stop button, overcharge management, IPX6 charger

Key Features

- Powered drive assistance for effortless patient transport
- Enhances workplace safety by reducing manual handling
- Intuitive force-sensing handles for proportional propulsion
- Retractable drive mechanism for manual push option
- Simple user interface with battery indicator, service alert & emergency stop
- IPX6 water resistant battery charger
- Designed exclusively for trained healthcare staff in clinical facilities

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au