

ALLEGRO WEIGHING DEVICE

SKU: WD03

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Allegro Weighing Device is a highly accurate and versatile tool designed to meet the diverse weighing needs in various settings, such as healthcare, fitness, and industrial environments. Known for its precision and reliability, this device plays a crucial role in applications where accurate weight measurement is essential.

FEATURES

- **High Precision:** The Allegro Weighing Device is engineered to provide highly accurate weight measurements, crucial for applications in healthcare for medication dosing and patient monitoring.
- **User-Friendly Interface:** Equipped with an intuitive interface, the device is easy to use, ensuring efficient operation even for those not familiar with technical equipment.
- **Durable Construction:** Built to withstand regular use in professional environments, the device is made from robust materials that guarantee longevity and consistent performance.
- **Versatile Functionality:** Ideal for a variety of settings, the device can be used for patient weighing in healthcare, measuring goods in industrial settings, or tracking weight in fitness and wellness programs.
- **Compact and Portable:** Designed to be compact, it can be easily transported and stored, making it suitable for facilities with limited space.
- **Digital Display:** Features a clear digital display that allows for easy reading of measurements, minimizing errors.
- **Battery Operated:** Typically powered by batteries, providing the convenience of portability and eliminating the need for constant power supply.

SPECIFICATIONS

- **Capacity:** 300kg
- **Graduations:** 100g
- **Dimensions:** 95mm (w) x 170mm (h) x 80mm (d)
- **Power supply:** Non-Rechargeable battery
- **Weight of scale:** 1kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au