ALRICK 2300 TILT CARE BED



SKU: 7-23T

PRODUCT DESCRIPTION

PRODUCT DESCRIPTION

The 2300T Series bed is the most versatile Hi-Low and tilt bed available for the healthcare industry. Featuring up to 300kg SWL sturdy frame, therefore being suitable for patients with a range of weight requirements, this bed offers sizes from single up to king widths.

The 2300T bed is excellent for in-home and facility care, as it has many options available, making this bed aesthetically pleasing and practical for both the patient and carers.

FEATURES

- Electric Hi-Low
- Manual Vascular Support
- Back lift handset
- APCO Grey & Ironstone powder-coated frame
- Curved natural oak head & footboards
- Electric kneebreak

- Electric backrest
- Forward and reverse Trendelenburg
- Bluetooth connectivity

*** Handset must be ordered separately ***

Standard inclusions

- Pull-out extension
- Auto chair position
- 125mm braking castors
- Underbed light

SPECIFICATIONS

2300 Tilt	Single	King Single	Extra Wide	Double	Queen	King
Overall Dimensions (mm)	2040 × 920	2040 × 1090	2040 × 1220	2040 × 1390	2040 × 1520	2040 × 1820
Sleeping Surface (mm)	2000 × 900	2000 × 1070	2000 × 1200	2000 × 1370	2000 × 1500	2000 × 1800
SWL (kg)	300	300	300	300	300	300
High/Low Position (mm)	830 / 390	830 / 390	830 / 390	830 / 390	830 / 390	830 / 390
Back Angle (degrees)	77	77	77	77	77	77
Knee Angle (degrees)	40	40	40	40	40	40
Reverse Angle (degrees)	16	16	16	16	16	16
Trendelenburg Angle (degrees)	16	16	16	16	16	16
Product Codes	23ST-000100	23KST-000100	23EWT-000100	23DT-000100	23QT-000100	23KT-000100

Brochure

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au