

AQUATEC OCEAN ERGO TRANSIT COMMODE

SKU: PKG-1643718-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Ocean Ergo Transit Commode, a thoughtfully designed aid for individuals with mobility challenges. This commode chair is engineered to provide a comfortable and safe bathroom experience, with features that cater to both the user's and caregiver's needs.

This kit also comes with a **Pan Holder** and **Soft Seat**.

FEATURES

- **Ergonomic Seat:** The ergonomically designed seat enhances comfort and accessibility, making personal hygiene tasks easier and more dignified.
- **Adjustable Footrests:** Footrests can be adjusted in height and angle to accommodate different user requirements, promoting proper posture and support.
- **Swing-Away Armrests:** For ease of transfer and to assist caregivers, the armrests swing away smoothly, providing unobstructed access.
- **Transit Wheels:** Equipped with four transit wheels, the commode can be easily maneuvered over toilets and into showers, with locks to ensure stability during use.
- **Corrosion-Resistant Frame:** Constructed with high-quality, corrosion-resistant materials suitable for use in wet environments.
- **Safety Features:** Includes brakes on all wheels and a stable base to ensure user safety while sitting or transferring.
- Tool-less height adjustable seat
- Machine washable backrest
- Space saving push bar
- Fold-away texture footplates

- Available in transit or self-propelled (only if you actually stock or can order both)

The Ocean Ergo Transit Commode stands out for its user-centric design, providing a safe and convenient solution for those who need assistance with daily toileting and showering tasks. Its blend of comfort, safety features, and ease of use makes it a valuable addition to any home or healthcare facility.

SPECIFICATIONS

Overall Width	565 mm
Overall Depth	905 mm
Seat Width	480 mm
Seat Depth	480 mm
Min Seat Height	475 mm
Max Seat Height	600 mm
Product Weight	17 kg
Max User Weight	Regular: 150 kg | Bariatric: 180kg
Brand	KCare

View the full Accessories Brochure

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au