

ARJO ALPHA ACTIVE 4 ALTERNATING PRESSURE SYSTEM

SKU: 648326-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Designed for the prevention and management of pressure injuries, the Alpha Active 4 pressure redistributing mattress replacement system is suitable for a range of care environments, including home and community care settings.

The Alpha Active 4 system has been designed as a cost-effective solution to aid in the prevention and management of pressure injuries.

The system offers advanced features for the delivery of high quality care to vulnerable patients and residents, whilst remaining easy-to-use for the carer and simple to clean and maintain.

The Arjo Alpha Active 4 system requires both the mattress replacement system and pump to be purchased.

FEATURES

- Electronic pressure setting allows accurate and appropriate pressure settings to be set for various weight ranges
- Variable cycle times
- AutoFirm provides a stable support surface for patient transfer and nursing procedures
- Two modes of operation: active / alternating and reactive / static
- Includes low pressure alarm, power fail alarm, and pump fault alarm
- Modular construction contributes to ease of cleaning and maintenance

- Ventilated standard air cells designed to aid airflow within the mattress
- Cell-on-Cell designed to provide additional support and to help increase patient stability
- Three static head cells
- Reliant IS2 cover is vapour permeable and water resistant, incorporating 2-way stretch and bacteriostatic properties
- Rapid deflation for CPR
- **Pump sold separately**
- **Both the Mattress Replacement System and Pump is required to operate the Arjo Alpha 4 Alternating Pressure System**

SPECIFICATIONS

Total Length	Total Width	Total Height	SWL
1950mm	850mm	200mm	200kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au