

ARJO BARIATRIC BASIC SLING

SKU: MAA8010-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Arjo Bariatric Basic Sling provides reinforced material and stitching in the most stressed areas of the sling, including the shoulders, legs and spine. The Bariatric Basic Sling has the same shape as the Deluxe Sling, but a more economical design. Padded leg sections provide extra support and assist with easy fitting and removal. Suitable for all types of general transfers.

FEATURES

- It provides reinforced material and stitching in the most stressed areas of the sling, including the shoulders, legs and spine.
- More Economical design than the deluxe sling.
- Suitable for all types of general transfers.
- Available in different sizes: M, L, XL, XXL

SPECIFICATIONS

- Safe working load: 454 kg/1000 lbs

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au