

# ARJO BATHING MESH SLING

**SKU:** MLA4060-H

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Arjo Bathing Mesh Sling** is designed to streamline patient bathing, showering, and personal hygiene routines. Its specialized mesh construction allows water to flow through freely and drain rapidly, making it ideal for wet environments and thorough washing. Featuring an unpadded leg design, this sling simplifies the process of dressing and clothing removal for hygiene routines, providing a smoother, more comfortable transfer experience for patients with spinal or skeletal conditions. With extended leg support for improved lower-body access and an integrated padded head support to assist individuals with reduced neck or trunk stability, it offers caregivers a safe, ergonomic, and dependable transfer solution.

### Features

- **Water-Permeable Mesh Fabric:** Enables rapid water drainage and fast drying, making it ideal for showering, bathing, and hygiene tasks
- **Unpadded Leg Sections:** Simplifies smooth application and effortless clothing removal during toileting and personal care
- **Integrated Padded Head Support:** Delivers essential comfort and upper-body stability for individuals unable to maintain head control independently
- **Extended Leg Reach:** Wraps securely around the thighs while providing caregivers with enhanced access to the lower body
- **High Weight Capacity:** Supports a Safe Working Load (SWL) of 272 kg for secure and dependable transfers
- **Broad Compatibility:** Fully compatible with 2-point and 4-point Arjo spreader bars and mobile/ceiling lifts

- **Durable Construction:** Engineered with robust, high-strength materials built to endure frequent laundering and rigorous daily care routines

## Specifications

| | |
|---|---|
| **Brand** | Arjo |
| **Product Name** | Bathing Mesh Sling (with Head Support) |
| **Base Article Number** | MLA4060 |
| **Available Sizes** | S, M, L, XL |
| **Product Codes** | MLA4060-S (Small)<br>MLA4060-M (Medium)<br>MLA4060-L (Large)<br>MLA4060-XL (Extra Large) |
| **Safe Working Load (SWL)** | 272 kg |
| **Material** | Heavy-duty, fast-draining mesh |
| **Design Configuration** | Extended leg loops (unpadded) with integrated padded head support |
| **Compatible Lifts & Spreader Bars** | Maxi Move: 2-point, 4-point<br>Maxi Twin: 2-point, 4-point<br>Maxi 500: 2-point<br>Minstrel: 4-point<br>Maxi Sky 440: Fixed 2-point<br>Maxi Sky 600: 2-point<br>Maxi Sky 2: 2-point, 4-point |
| **Application** | Patient transfers during bathing, showering, and toileting routines |

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

**Brighton Showroom**
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

**Elizabeth Vale Showroom**
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

**Reynella Showroom**
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

## VIC

**Frankston Showroom**
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

**Monash Showroom**
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

## NT

**Darwin Showroom**
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au