

ARJO CAREVO - HEIGHT ADJUSTABLE SHOWER TROLLEY

SKU: BAC1113-01

PRODUCT DESCRIPTION

ARJO CAREVO PURPOSE & BENEFITS

The Arjo Carevo shower trolley is designed to provide a comfortable shower experience for patients and residents who need to remain in a lying position, while enabling the caregiver to work efficiently and ergonomically.

Carevo is height-adjustable and features patented Ergo-Access areas to enable the caregiver to get closer to the resident and maintain a comfortable working posture, while the patented Flexi Zone is designed to enhance comfort for the patient or resident.

FEATURES

- Height adjustable shower trolley
- Facilitates improved working posture
- Clearly marked control panels located on both sides and at the head of the trolley, allowing easy access for the carer
- The soft surface material and patented Flexi Zone are designed to enhance patient or resident comfort throughout the hygiene routine
- Comfort handles provide support for the resident or patient
- Purchase includes: head pillow | battery (x2) | battery charger | battery holder | drainage hose – upper drainage | lower drainage
- Hire includes: head pillow | one battery | battery charger | battery holder | drainage hose – upper drainage | lower drainage

SPECIFICATIONS

Total Length	2005mm
Width | Stretcher	586mm
Total Width | Including Side Supports	759mm – 888mm
Chassis Length	1211mm
Chassis Width	700mm
Height | Top Position Floor to Stretcher	955mm
Height | Low Position Floor to Stretcher	605mm
Stroke Range	350mm
Weight | Including Mattress	82kg
Maximum Patient Weight	182kg
SWL | Patient and Water	200kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au