

# ARJO GENERAL PURPOSE POLYESTER SLING WITH HEAD SUPPORT

**SKU:** MLA2000-H

## PRODUCT DESCRIPTION

## PURPOSE & BENEFITS

Arjo standard padded loop sling is used on 2-point loop spreader bars. The slings are intended for everyday patient handling routines. The padded sling with head support offers padding for the leg area of the sling, increasing comfort and helping to prevent the material from roping and creasing under the resident's/patient’s legs. Padded loop slings are available in four sizes (S, M, L and XL).

## FEATURES

- It provides good comfort and support for everyday residents/patients handling routines, such as bed and chair-related transfers.
- Padded sling with head support.
- Increased comfort and helps to prevent the material from roping and creasing under the resident's/patient’s legs.
- Available in different sizes: S, M, L, XL

## SPECIFICATIONS

- Safe working load: 272 kg/600 lbs

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

**Brighton Showroom**
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

**Elizabeth Vale Showroom**
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

**Reynella Showroom**
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

## VIC

**Frankston Showroom**
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

**Monash Showroom**
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

## NT

**Darwin Showroom**
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au