

# ARJO MINUET

**SKU:** 100500AUS

## PRODUCT DESCRIPTION

## PURPOSE AND BENEFITS

Designed for use in residential and convalescent homes and private residences, our range of Minuet 2 community beds provides a real solution for independent living. Advanced electrical profiling, variable bed height, combined with a domestic appearance makes the Minuet 2 the natural choice for the private, residential and convalescent home user. To facilitate transport and storage the bed has been designed to be demounted and foldable to minimise space requirements.

**View the Arjo Minuet 2 product brochure**

## FEATURES

- Low platform height 250 mm (10”) from the floor
- Safe working load 180 kg (396 lbs)
- Electric backrest, knee break and height adjustment
- Pro-Contour Advance profiling feature has demonstrated low interface pressures which may help in the prevention of pressure injuries
- Battery back-up.
- Integrated mattress platform extension of 130 mm (5”) with mattress retainer
- Designed for use in private, residential and convalescent homes
- Easy to disassemble into four parts for improved installation access
- Optional two-way tilt; foot-down/head-down*
- Compact storage/transportation (using transit brackets)
- Wide range of accessories available.

## SPECIFICATIONS

| **160 & 161** | |
|---|---|
| Overall length | 2250 mm (88.6”) |
| Overall length (extended) | 2390 mm (94”) |
| Overall width without safety sides | 940 mm (37”) |
| Overall width with safety sides | 1050 mm (41.3”) |
| Mattress size (L x W x H) | 1980 x 860 x 125 mm (78 x 33.8 x 5”) |
| Height range | 250 – 660 mm (9.8 – 26”) |
| Tilt range (if fitted) | +/– 10.8° |
| Backrest angle (maximum) | 69° |
| Thigh elevation (maximum) | 32º |
| Safe working load | 180 kg (396 lbs) |
| Maximum patient weight | 145 kg (319 lbs) |
| Section weights - Head end | 22.8 kg (50 lbs) |
| Section weights - Foot end | 22.8 kg (50 lbs) |
| Foot end mattress platform section | 19.8 kg (44 lbs) |
| Head end mattress platform section | 23.5 kg (52 lbs) |
| Electrical data (all models) | |
| Power in | 1.7 A max. 120 V a.c. 60 Hz / 1 A max. 230 V a.c. 50 Hz |
| Electric shock protection | Class II, Type B |
| Ingress protection | IP24 |
| Battery | 12 V, 1.2 Ah sealed rechargeable lead/acid gel |

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

**Brighton Showroom**
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

**Elizabeth Vale Showroom**
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

**Reynella Showroom**
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

## VIC

**Frankston Showroom**
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

**Monash Showroom**
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

## NT

**Darwin Showroom**
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au