

ASPIRE BED SELF HELP POLE

SKU: BEA013100

PRODUCT DESCRIPTION

FEATURES

- Free standing support handle provides assistance for the user to reposition or transfer in and out of bed
- Handle height is easily adjusted to suit the user
- Sits behind or beside the bed
- Leg length: 750mm
- Leg width: 800mm
- Leg height (clearance needed): 30mm
- Handle height minimum: 910mm
- Handle height maximum: 1450mm
- Safe working load: 150kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au