

## ASPIRE SHELL CHAIR

**SKU:** CHP19891-H

## PRODUCT DESCRIPTION

## PURPOSE & BENEFITS

The Aspire Shell Chair is a comfortable mobile seating solution designed to maximise user access and social interaction within an Aged Care Facility or Healthcare environment. The compact design allows access under tables for sharing meals, games and other activities. The lightweight frame with anterior tilt function for sit-to-stand transfers reduces the manual handling burden for caregivers.

## FEATURES

- Forward tilt functionality assists with chair egress and standing transfers
- Backwards tilt functionality assists in pressure redistribution and patient stability during transfers
- Comfort & Support - Ergonomic wing headrest and lumbar cushioning for improved seating posture
- Storage Pocket - Conveniently carries everyday items such as clipboards, magazines or diaries
- Lightweight and Manoeuvrable - 12.5” rear wheels with attendant foot brake, 5” front swivel castors together with a lightweight design, makes for easy and safe handling for caregivers and a smooth and comfortable ride for users

## SPECIFICATIONS

| Display Name | Seat Depth | Seat Width | Seat To Floor Height | Safe Working Load | Unit Weight |
|---|---|---|---|---|---|
| Aspire Shell Chair - 450mm Wide | 440 mm | 450 mm | 450 - 530 mm | 135 kg | 30 kg |

| | | | | | |
|---|---|---|---|---|---|
| Aspire Shell Chair - 510mm Wide | 470 mm | 510 mm | 450 - 530 mm | 135 kg | 30 kg |

## BROCHURE

Aspire Shell Chair Brochure

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

**Brighton Showroom**
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

**Elizabeth Vale Showroom**
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

**Reynella Showroom**
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

## VIC

**Frankston Showroom**
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

**Monash Showroom**
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

## NT

**Darwin Showroom**
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au