

ASPIRE VOGUE FOREARM SEAT WALKER

SKU: WAF705400

PRODUCT DESCRIPTION

Purpose & Benefits

The **Aspire Vogue Forearm Seat Walker** is designed for users who require additional upper-body support while walking. Featuring contoured forearm pads and advanced ergonomic adjustments, it promotes correct posture, balance, and comfort while providing smooth and stable mobility indoors or outdoors.

Features

- Soft, contoured forearm pads ensure maximum comfort and support
- Smart Posture System allows height, depth, and angle adjustment for optimal balance and weight distribution
- Danish-inspired design with ergonomic handles, profiled tubing, and refined colour finish
- Lightweight frame and small turning circle for superior manoeuvrability
- “AutoLocate” system enables smooth and precise height adjustment to pre-set positions
- Locking mechanism secures the frame when folded for easy transport and storage
- Advanced locking handbrakes provide additional safety when walking or seated

- Unbreakable TPE super soft wheels ensure durability and a smooth, comfortable ride

Specifications

- **Seat Width:** 441 mm
- **Total Height:** 900 – 1170 mm (adjustable)
- **Seat to Floor Height:** 620 mm
- **Safe Working Load (SWL):** 150 kg
- **Unit Weight:** 10 kg (including accessories, bag, and backstrap)
- **Frame Weight:** 9.1 kg
- **Colour:** Black
- **Warranty:** 1 year

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au