

ASPIRE VOGUE INDOOR WALKER WITH TRAY

SKU: WAF705100-H

PRODUCT DESCRIPTION

Purpose & Benefits of the Aspire Vogue Indoor Walker With Tray

The Aspire Vogue Indoor Walker is designed to provide mobility support for individuals who require assistance with walking in indoor environments. This stylish and functional walker is ideal for use at home or in residential care settings, offering features that enhance safety, convenience, and independence.

Features

- **Sleek Design**: Features a modern and stylish design that fits seamlessly into home environments, providing a more attractive alternative to traditional walkers.
- **Compact and Manoeuvrable**: Designed with a compact frame and swivel front wheels, making it easy to navigate through tight spaces, narrow hallways, and around furniture.
- **Convenient Storage Tray**: Includes a built-in storage tray that allows users to carry personal items, meals, or other essentials easily, enhancing independence and convenience.
- **Easy Folding Mechanism**: Features an easy-to-use folding mechanism that allows the walker to be quickly folded for storage or transport, making it ideal for use in small living spaces or when travelling.
- **Ergonomic Hand Grips**: Comes with comfortable, ergonomic hand grips that provide a secure and comfortable hold, reducing hand fatigue during use.
- **Stable and Durable**: Constructed from high-quality materials to ensure stability and durability, providing reliable support for everyday use.

Specifications

- **Safe Working Load**: 110 kg

- **Dimensions**:
  - Height: 860–970 mm
  - Length: 650 mm
  - Width: 540 mm
- **Colour Options**: Red, Champagne
- **Weight**: 5.1 kg

The Aspire Vogue Indoor Walker is an excellent choice for individuals seeking a stylish, practical, and reliable mobility aid for indoor use. Its combination of modern design, convenient features, and user-friendly functionality makes it a valuable tool for enhancing mobility and independence in the home.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au