

BETTERLIVING COMMUNITY BED

SKU: BE1005-H

PRODUCT DESCRIPTION

Purpose & Benefits

The ideal choice for community or aged care settings, offering a range of optional accessories, ease of assembly, and straightforward storage, making it a versatile and adaptable solution. The BetterLiving Community Bed functions as a four-part profiling bed. This wide range of adjustments ensures that the bed can be positioned in multiple ways to cater to various needs and preferences. The bed incorporates a built-in base extension that increases the length of the bed platform, providing added comfort. It comes equipped with a lightweight and water-resistant hand controller, featuring soft-touch buttons, including a one-touch auto-contouring function.

The autoregressive backrest helps minimise pinching or shear. The BetterLiving Community Bed is powered by a high-quality and reliable Linak electric drive system. The mattress platform is easy to clean and service, and it includes removable mattress retainer handles for added convenience. The bed is equipped with four high-quality individual locking castors that enable easy manoeuvrability.

Features

- Electric hi-lo height adjustment, adjustable backrest, knee break, and Trendelenburg tilt, allowing for multiple positioning options
- Lightweight and water-resistant hand controller with soft touch buttons
- One-touch auto-contour function is ideal for patients with cognitive impairment
- Auto contour function simultaneously raises both the head and knees to ensure optimal positioning, which reduces shear when the backrest is raised and the mattress folds
- Powered by a high-quality Linak electric drive system for durability and dependable performance

- Easy to clean mattress platform
- Features a built-in extension, increasing the bed length by 200 mm for added comfort
- Four individually locking castors

Bed Sizes

- Single Bed
- King Single Bed

Optional Accessories

- BetterLiving Community Bed, Side Rail – Three Quarter Length (Pair)
- BetterLiving Self Help Pole
- BetterLiving Grab Bar
- BetterLiving IV Pole

Specifications

	Single	King Single
Total Length	2170mm	2170mm
Total Length | Extended	2370mm	2370mm
Total Width	995mm	1220mm
Recommended Mattress Size	1960mm x 865mm x 170mm max.	1960mm x 1050mm x 170mm max.
Height	330mm – 730mm	330mm – 730mm
Maximum Patient Weight	180kg	180kg
SWL	220kg	220kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au