

BETTERLIVING ELEPHANT FEET BED RAISERS & CHAIR RAISERS - SET OF 4

SKU: H-BL16

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

These plastic Elephant Feet Bed Raisers & Chair Raisers are designed to go under the legs of your bed, lounge chair, couch or table to increase the overall height. They allow you to ensure the height of your furniture is right for your needs. Bed Raisers can also be used to raise one end of your bed, to create a slight incline, which can be helpful in managing some digestive, respiratory and circulatory problems.

Made from durable, high impact plastic, these stackable Bed Raisers are available in three different heights. The plastic raiser has a top well that can accommodate square or round legs up to 65 mm wide. They are sold in sets of four. Each set of four can hold up to 250 kilograms of weight, including the furniture and any people.

Please note that your furniture must be sufficiently heavy to remain stable on top of the raisers, so these products are not suitable for lightweight tables or chairs.

The 40mm high raiser has a base width of 125mm. The 100mm high raiser has a base width of 198mm. Combine these two to create a 140mm high raiser with a base width of 198mm.

KEY FEATURES

- Durable plastic furniture raisers
- Place under furniture legs to increase the overall height
- Use for beds, couches, tables or chairs
- Suits square and round legs
- Sold in sets of 4

- Holds up to 250kg total

You can buy these elephant legs for beds online using AfterPay and NDIS payment methods.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom
343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom
Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112
08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom
Shop 1/141 Sherriffs Rd, Reynella SA 5161
08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom
111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451
frankston@medimart.com.au

Monash Showroom
23 Winterton Road, Clayton VIC 3168
03 7503 8354
monash@medimart.com.au

NT

Darwin Showroom
Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468
support@medimarnt.com.au